

**SEE BETTER
PERFORM
BETTER**
Dental Products

Choosing the right frame for your loupes is just as important as choosing the right magnification!



What Sets Enova Loupes Apart

- 1 Lightweight, aircraft aluminum construction
- 2 Adjustable silicone nose pads for added comfort; available in three styles
- 3 Magnetic light adapter for quick and easy headlight attachment
- 4 TriVEX® lenses for durability and optical performance
- 5 Spring-loaded temple arms to maintain telescope convergence and a secure, comfortable fit

Fast fact: As magnification increases, the field of view decreases. Many first-time loupe buyers prefer 2.5x magnification to give them a larger depth of field and field of view.

Loupes

Vizix® Loupes

Galilean | 2.5x, 3.0x, 3.5x

Colors: Silver, Green, Black, Purple, and Blue



Prismatic | 3.5x, 4.5x, 5.5x

Color: Black/Silver



Ergo Elite | 3.5x, 4.5x, 5.5x

Color: Black/Silver



Frames

Durable Frame Options Suited to Your Style

Enova knows that quality, durability, fit, and aesthetics are important aspects when it comes to selecting a frame that works with your loupe and lighting solution. We have taken this into account with a curated collection of aluminum frames that work seamlessly with our loupe and lighting products.

Airon® Frame, Regular or Small

Colors: Blue, Black, Silver, Red, and Pink



Frames

Airon® Radius Frame

Colors: Blue, Black, Silver, Red, and Pink



Flak® 2.0 XL Frame

Colors: Black, White



Centerboard® Frame

Color: Clear



Note: Some frames will not accommodate high power prescriptions.

Loupe Lights



AXIS

Lightweight Light

Illumination: 30,000 LUX

Color Temperature: 6,100K

Spot Size: 4" (10.2cm) at 12" (30cm) distance

Control: 3 level brightness control

Total Weight: 1.1oz (30g)

Run-Time: 10+ hrs with a PP1.5



ATOM PLUS

Cordless Light

Illumination: 27,000 LUX

Color Temperature: 6,000K

Spot Size: 3" (7.6cm) at 12" (30cm) distance

Touch Sensor: 3 level brightness control,
power on/off

Total Weight: 1.4oz (41g)

Run-Time: 3 hrs

Loupe Light Accessories

Activ™ Neck Battery provides:

- 5+ hours of working time
- Capacitive touch technology with three brightness settings
- Color options including black, purple, blue, and green
- Compact, Lightweight Design:
1"x1" (2.5cm x 2.5cm)
5oz (149g)



PP1.5 Battery provides:

- Rechargeable Lithium-ion battery
- 10+ hours of working time
- Battery pack hip belt
- Compact, Lightweight Design:
4"x2.5"x1" (10.2cm x 6.6cm x 2.5cm)
7oz (200g)



Fast fact: Using a headlight to illuminate the working field not only improves procedure performance, but can also lengthen your career. LED lighting, coupled with magnification, improves ergonomic posture and reduces musculoskeletal disorders that often worsen over time.





ENOVA SATISFACTION GUARANTEE

Use our products and love them. If for any reason, you're not satisfied with our loupes or lighting, contact your Aluro representative for a visit. Go to our website to view Terms & Conditions

To learn more about how Enova can help you see better, so you can perform better...

Visit www.aluro.co.nz / Call 0800 800 085 / Email orders@aluro.co.nz

Healthcare
Aluro[®]
NZ
making dentists smile